

A2J GmbH

Technology that connects

For all uses
WATERPROOF



**Stability and comfort
in 3D - technology**

PELVI.LOC®



Technology that connects

PL-3DWP-BL-L

Waterproof pelvis retraction belt
in black



Design

Special closed-pore material, impermeable to water and water repellent! The pelvic pads are 3-dimensional for the optimum positioning of the pelvis without pressure on the anterior superior spine. The fit of the belt is set with easy-to-use ratchet buckles made from stainless steel. The punched straps on the frame are length-adjustable by up to 8cm, which means they ensure a situation-adapted, room for manoeuvre. For the securing of the wheelchair user, the pelvic pads are connected in the middle with a toothed strap.

Item no.

PL-3DWR-BL-L for SB 36-46 cm

Use

For the positioning of the pelvis in the shower wheelchair, bathingsupport, active wheelchair, rehab buggy, seat shells and much more.

Special attributes

Faceted thinning out of the padding material in the transition from the pad to the fastening belt. Prevents lateral pressure points on the pelvis. V-R-24-BL-ERGO-170 toothed strap

PL-3DWP-WH-L

Waterproof pelvis retraction belt
in white



WATERPROOF

Additional item:

V-R-24-BL-ERGO-225 V toothed strap, long and riveted
V-R24-BL-LR V punched strap, long and riveted

Tube brackets for attaching the belt to the frame of the wheelchair are available in the following diameters:

Item no.

PL-3DWR-WH-L for SB 36-46 cm

Item no.

RS-20	tube diameter 20mm
RS-22	tube diameter 22mm
RS-25	tube diameter 25mm
RS-30	tube diameter 30mm
RS-32	tube diameter 32mm

A2J GmbH

Technology that connects



WATERPROOF

Subject to technical modifications and errors. Variations in colour from the depicted versions are possible.