
		Price 2019 excl. 19% VAT Model: Crazy Horse Active	
		<input type="checkbox"/> Offer <input type="checkbox"/> Order	
		Date:	
		Name:	
		Commission:	
		Billing address:	
		Company:	
		Name:	
		Street:	
		City:	
		Post code:	
Phone:			
Email:			
Customer details:		Delivery address:	
Wheelchair:		Company:	
Seat width:	cm	Name:	
Seat depth:	cm	Street:	
Lower leg length:	cm	City:	
Manufacturer:	Model:	Post code:	
Height:	Body weight:	Phone:	
Thru axle diameter:	<input type="checkbox"/> 12,00mm <input type="checkbox"/> 12,7mm	Email:	
Base price:			
Crazy Horse Active			755,46 €
Art.-no.:	Color	Weight in grams	Price
700001	Matt black anodized aluminum	-	free of charge
700002	Matt silver anodized aluminum	-	free of charge
Rim (spoke + hub + brake disc)			
720001	30 mm wide matt silver	1020	default
Tires (tire + tube + flap)			
710001	Innova 20 x 2.00 black city	860	default
710002	20 x 2,125 white wall city	720	32 €
710003	Schwalbe 20 x 2,250 jumpin jack black MTB	860	34 €
710004	SPO 20 x 3.00 black sand	1620	38 €
710005	SPO 20 x 3.00 white wall sand	1620	39 €
710006	SPO 20 x 3.45 black sand	1400	49 €
710007	SPO 20 x 4 ¼ black sand	2360	49 €
Handlebar			
740001	Handlebar individually adjustable in angle	-	35 €
Parking brake			
740101	Brake lever Magura lockable	-	29 €
Fender			
740601	80mm fender chromed	-	48 €
Subtotal net excl. VAT			

ACCESSORIES:

Art.-no.:	Change wheel rim / 64 mm (w/o tupe and tire 1680 gr.)	Weight in grams	Surcharge
720101	64 mm rim matt black Innova 20 x 2.00 city	2520	148 €
720102	64 mm rim matt black Schwalbe 20 x 2.25 MTB	2540	154 €
720103	64 mm rim matt black Kenda 20 x 2.175 MTB	2380	149 €
720104	64 mm rim matt black Kenda 20 x 2.125 city white wall	2380	152 €
720105	64 mm rim matt black GGG 20 x 3.45 sand	3100	169 €
720106	64 mm rim matt black SPO 20 x 3.00 sand	3300	148 €
720107	64 mm rim matt black SPO 20 x 3.00 white wall sand	3300	149 €
720108	64 mm rim matt black SPO 20 x 4 ¼ black sand	4000	169 €
	Change wheel rim / 80 mm (w/o tupe and tire 1940 gr.)		
720201	80 mm rim matt black Innova 20 x 3.00 black sand	3540	168 €
720202	80 mm rim matt black Innova 20 x 3.00 white wall sand	3540	169 €
720203	80 mm rim matt black GGG 20 x 3,45 black sand	3540	179 €
720204	80 mm rim matt black SPO 20 x 4 ¼ black sand	4000	179 €
	Wheels with axle		
772112	MTB colour black 25 Zoll axis 12,0 mm x 104 mm		450 €
772111	MTB colour black 25 Zoll axis 12,7 mm x 104 mm		450 €
772109	Flame colour white wall 25 Zoll rim 64 mm axis 12,0 mm x 104 mm		750 €
772110	Flame colour white wall 25 Zoll rim 64 mm axis 12,7 mm x 104 mm		750 €
772114	Flame colour white wall 25 Zoll 80 mm axis 12,0 mm x 104 mm		950 €
772113	Flame colour white wall 25 Zoll 80 mm axis 12,7 mm x 104 mm		950 €
772115	Hot Frog III colour black 26 Zoll wheel 80mm hole-rim 12,0 mm x 104 mm		950 €
772116	Hot Frog III colour black 26 Zoll wheel 80mm hole-rim 12,7 mm x 104 mm		950 €
	Flap for rim with holes		
710300	Flap 80 mm wide matt black		free of charge
710301	Flap 80 mm wide matt white		free of charge
710302	Flap 80 mm wide matt orange		free of charge
710303	Flap 80 mm wide matt blue		free of charge
710304	Flap 80 mm wide matt silver		free of charge
	Puncture protection from Teflon rear wheels		
720113	1 pair for insertion between tube and tire		28 €
	Rearview mirror		
740700	rearview mirror left		35 €
740701	rearview mirror right		35 €
	Rear light		
740201	rear light		88 €
	Air pump with battery		
740204	High pressure pump with battery		78 €
740205	Car charger for high pressure pump		16 €
740206	0.5m extension hose with clip		12 €
	Freight / Assembly		
600001	Shipping on pallet with insurance (Ask for daily rates)		on request
600002	Factory pickup A2J GmbH		free of charge
600003	Assembly of the wheelchair Crazy HorsE-Power in the factory A2J GmbH		95 €



Technik die verbindet

Price:	
Subtotal:	
+19% VAT	
Total price incl. 19% VAT:	

Comments:

Privacy Regulation:
<input type="checkbox"/> I hereby consent to the processing of my personal data in accordance with DSGVO. See privacy policy A2J website.

Legal notice at: 6km / h - 15km / h -30km / h - open power Crazy-HorsE-Power
<input type="checkbox"/> I have read and understood the legal notice on the A2J website.

The prices are exclusive of the statutory 19% VAT.

* The Crazy HorsE-Power is a medical device according to MDD 93/42 / EWG.

** Important note: Use only on private property, abroad please inform yourself about local regulations.