

PELVI.LOC® PL-2DA Active

PELVI.LOC® PL-2D Thigh Belts

HMV.Nr.: 18.99.99.0905



Version

2D-belt for enhanced sports activities with enlarged pads - surface for improved pressure distribution.

The padded covers give maximum protection of participants and prevent buckle from releasing in use. Hole strap on frame adjustable up to 6 cm.

Special features

Monoski racer with Martin Braxenthaler:

- **PL-2DA-M** as a chest strap and thigh strap



Amendment

Tube clamps for attaching the belt to the wheelchair frame available in the following diameters:

Article numbers

RS-20: Tube diameter 20 mm

RS-22: Tube diameter 22 mm

RS-25: Tube diameter 25 mm

RS-30: Tube diameter 30 mm

RS-32: Tube diameter 32 mm



Scope of application

- Safe driving with 2D pelvic positioning in a wheelchair.
- Can be used for thorax positioning in a standing device.
- Pelvic positioning with pressure distribution and stability in sports wheelchairs.

Special features

Article numbers




PL-2DA-S for SW 20-32 cm

PL-2DA-M for SW 30-40 cm

PL-2DA-L for SW 36-46 cm

Technical Data

	XXXX Bambino	XXS	XS	S	M	L	XL
Seat width (cm)				20 - 32	30 - 40	36 - 46	
Padding size (cm)				19 x 11	23 x 12	28 x 14	
Weight (g)				220	280	280	
Predetermined breaking point (at approx. kg)				60	80	80	

Ratchet size XXXS - XXS			
Ratchet size XS - S	2,4 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-24-FLAT-LO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap	LR-BL-320	Produkt ist latexfrei 	LIZENZIERTER PARTNER REG. NR. DE-59142/47496739-V

PELVI.LOC® PL-2DAS Active for moulded Seats

PELVI.LOC® PL-2D Thigh Belts

HMV.Nr.: 26.99.02.2000



Version

Thigh and pelvis-retraction belt also for moulded seats and Active- Wheelchairs.

Maximum coverage with contoured form for enhanced sports activities. The padded covers give maximum protection of participants and prevent ratchet from releasing in use.

Performance


With other Buckles (B-24-Star) available.
Picture: EM Frankfurt 2013, Wheelchair basketball, national team



Amendment

Set of securing buckles for molded seat belt for attachment to the seat shell or wheelchair.

Article number
PB-BL-53



Scope of application

For all wheelchair sports activities, where secure fastening to the chair is an advantage.

For example: Wheelchair Basketball, Wheelchair Rugby, Tennis, Table tennis, Electric Wheelchair hockey, Ice hockey, hand cycling, and not to forget winter sports i.e. mono ski and cross country skiing.

Special features

Formed shape with ultra thin end at the bucket seat fixing point to avoid lateral pressure marks at the pelvis.

Article numbers

PL-2DAS-S for SW 20-32 cm

PL-2DAS-M for SW 30-40 cm

PL-2DAS-L for SW 36-46 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)				20 - 32	30 - 40	36 - 46	
Padding size (cm)				19 x 11	23 x 12	28 x 14	
Weight (g)				220	280	280	
Predetermined breaking point (at approx. kg)				60	80	80	

Ratchet size XXXS - XXS		<p>Product is latex-free</p>	
Ratchet size XS - S	2,4 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-24-FLAT-LO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap	LR-BL-320		

PELVI.LOC® PL-2DB

PELVI.LOC® PL-2D Thigh Belts

HMV.Nr.: 18.99.99.0905



Version

2DB belt suitable for thigh. The buckles are easily accessible

The ratchet will run in addition to a strap to prevent injuries. Hole strap on frame adjustable up to 6 cm. For fixing on to the frame use compatible clamping clips.

Performance

Cross country skis with:

- **PL-3DSRR-M** 3D moulded seat Reverso with twisted toothed belt and belt set
- **PL-2DB-XS** Thigh Belt Front



Amendment

Tube clamps for attaching the belt to the wheelchair frame available in the following diameters:

Article numbers

- RS-20:** Tube diameter 20 mm
- RS-22:** Tube diameter 22 mm
- RS-25:** Tube diameter 25 mm
- RS-30:** Tube diameter 30 mm
- RS-32:** Tube diameter 32 mm



Scope of application

- For pelvis positioning in a wheelchair and sports wheelchair



Article numbers

- PL-2DB-XXS** for SW 12-20 cm
- PL-2DB-XS** for SW 16-24 cm
- PL-2DB-S** for SW 20-32 cm
- PL-2DB-M** for SW SB 30-40 cm
- PL-2DB-L** for SW SB 36-46 cm
- PL-2DB-XL** for SW 44-52 cm



Technical Data

	XXXS	XXS	XS	S	M	L	XL
Seat width (cm)		12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)		12 x 7	17 x 7	20 x 7	22 x 8	25,5 x 8	28,5 x 9
Weight (g)		200	220	220	280	280	280
Predetermined breaking point (at approx. kg)		40	60	60	80	80	8

Ratchet size XXXS - XXS	1,4 x 0,2 x 14,0 cm	Product is latex-free 	 LIKENZIERTER PARTNER REG. NR. DE 591427/4696739-V
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS	B-14-STAR-BL		
Buckle size XXS	B-14-GO-BL		
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		

PELVI.LOC® PL-2DS Basic for moulded Seats

PELVI.LOC® PL-2D Thigh Belts

HMV.Nr.: 18.99.99.0905



Version


2DB-belt suitable for thigh.

The buckles are easily accessible. The ratchet will run in addition to a strap to prevent injuries. Webbing in sizes XS, S, M and L

Amendment

For attachment to the seat shell or on the wheelchair.

Article number
PB-BL-38



Scope of application



- For pelvis positioning in a wheelchair and sports wheelchair.

Article numbers

PL-2DS-XS for width SB 16-24 cm
PL-2DS-S for width SB 20-32 cm
PL-2DS-M for width SB 30-40 cm
PL-2DS-L for width SB 36-46 cm

Technical Data

	XXXS	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	
Padding size (cm)			18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	
Weight (g)			220	220	280	280	
Predetermined breaking point (at approx. kg)			60	60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	 LIKENZIIERTER PARTNER REG. NR. DE 5914274696739-V
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		

PELVI.LOC® PL-3DA Activ

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

3D-active pelvis retraction belt

3D-belt high tensile forces with enlarged pads - top. Pelvic surface. 3-Dimensional shape gives maximum support of the pelvis without pressure on the spine "anterior superior". The covered buckles help protect from injury and accidental opening during use. Hole strap on frame adjustable up to 6 cm.

Amendment

Clamping clips for wheelchair frame. Other sizes available on request!

Article numbers

RS-20 Tube diameter 20 mm
RS-22 Tube diameter 22 mm
RS-25 Tube diameter 25 mm
RS-30 Tube diameter 30 mm
RS-32 Tube diameter 32 mm



Scope of application



Extensive 3D-belt for high tensile strength (e.g. spasmodic, wheelchair-sports). Covered buckle avoids the unintentional opening of the belt.

Article numbers

PL-3DA-XS for width 16-24 cm
PL-3DA-S for width 20-32 cm
PL-3DA-M for width 30-40 cm
PL-3DA-L for width 36-46 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	
Padding size (cm)			17,5 x 9	20 x 11	23 x 13	25,5 x 14	
Weight (g)			220	220	280	280	
Predetermined breaking point (at approx. kg)			60	60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	 LIZENZIIERTER PARTNER REG. NR. DE5914274696739-V
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			

PELVI.LOC® PL-3DAS Activ for moulded Seat

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Pelvis retraction belt for moulded seats

3D-belt high tensile forces with enlarged pads -surface. Nylon strap for passing through the moulded seat. 3-Dimensional shape gives maximum support of the pelvis without pressure on the spine "anterior superior". Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis. The covered buckles help protect from injury and accidental opening during use.

Amendment


For moulded seat for better pelvic positioning in the wheelchair.

Article number: PB-BL



Attachment kit for seat buckets.

Article number: BS-24-BL



Scope of application

For custom-made seating systems, moulded seat.

Special features

To improve positioning of the pelvis - retraction of the seat, with high spasticity.

Article numbers

PL-3DAS-XS for width 16-24 cm

PL-3DAS-S for width 20-32 cm



PL-3DAS-M for width 30-40 cm

PL-3DAS-L for width 36-46 cm

PL-3DAS-XL for width 44-52 cm

Technical Data

	XXXX Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)			17,5 x 9	20 x 11	23 x 13	28 x 14	29 x 16
Weight (g)			220	220	280	280	280
Predetermined breaking point (at approx. kg)			60	60	80	80	80

Ratchet size XXXS - XXS		<p>Product is latex-free</p> 	
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Ratchet size XL	2,4 x 0,2 x 22,5 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			

PELVI.LOC® PL-3DASR Activ reversed for moulded Seat

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Reversed pelvis retraction belt for moulded seat

Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis. 3-Dimensional shape gives maximum support of the pelvis without pressure on the spine "anterior superior".


Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis. The covered buckles help protect from injury and accidental opening during use.

Amendment

Pipe clamps for attaching the belt to the wheelchair frame available in the following diameters.

Article numbers

RS-20 Tube diameter 20 mm
RS-22 Tube diameter 22 mm
RS-25 Tube diameter 25 mm
RS-30 Tube diameter 30 mm
RS-32 Tube diameter 32 mm



Scope of application

For custom-made seating systems, moulded seat.

Special features



On the moulded seat remain attached the ratchets. Remove to improve positioning of the pelvis - retraction of the seat, with high spasticity.

Article numbers

PL-3DASR-XS for width 16-24 cm
PL-3DASR-S for width 20-32 cm
PL-3DASR-M for width 30-40 cm
PL-3DASR-L for width 36-46 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	
Padding size (cm)			17,5 x 9	20 x 11	23 x 13	25,5 x 14	
Weight (g)			220	220	280	280	
Predetermined breaking point (at approx. kg)			60	60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	
Ratchet size XS - S	1,9 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			

PELVI.LOC® PL-3DB Basic

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Pelvis retraction belt

The pelvic supports are 3-dimensional shape for optimal fixation of the pelvis without pressure on the spine "anterior superior ". The ratchet can also be attached with a tab to avoid injury. Hole strap on frame adjustable up to 6 cm.

Amendment

Tube clamps for attaching the belt to the wheelchair frame available in the following diameters.

Article numbers

RS-20 Tube diameter 20 mm
RS-22 Tube diameter 22 mm
RS-25 Tube diameter 25 mm
RS-30 Tube diameter 30 mm
RS-32 Tube diameter 32 mm



Scope of application



Wheelchair, reha-buggy, moulded seat.

Article numbers

PL-3DB-XS for width 16-24 cm
PL-3DB-S for width 20-32 cm
PL-3DB-M for width 30-40 cm
PL-3DB-L for width 36-46 cm
PL-3DB-XL for width 40-48 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)			18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)			220	220	280	280	280
Predetermined breaking point (at approx. kg)			60	60	80	80	80

Ratchet size XXXS - XXS		Product is latex-free 	 LIZENZIERTER PARTNER REG. NR. DE5914274696739-V
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			

PELVI.LOC® PL-3DS for moulded Seat

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

The pelvic pads are 3 - dimensional shape, for optimum fixation of the pelvis without pressure on the spine anterior superior. Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis.

Prevents lateral pressure points on the pelvis. Nylon strap for passing through the moulded seat. The ratchet will run in addition to a strap to prevent injuries.

Performance


Seat shell with **PL-3DS-S** pelvic retraction belt.



Amendment

For mounting on the moulded seat or on the wheelchair.

Article number: **PB-BL**



Scope of application

For custom made seats/ moulded seats.

Special features


Facet-like thinning of the padding material in the transition from the pad to the fastening strap. Prevents lateral pressure points on the pelvis.

Article numbers

PL-3DS-XXXX for SB 8-16 cm Bambino
PL-3DS-XXS for SB 12-20 cm
PL-3DS-XS for SB 16-24 cm
PL-3DS-S for SB 20-32 cm
PL-3DS-M for SB 30-40 cm
PL-3DS-L for SB 36-46 cm
PL-3DS-XL for SB 44-52 cm

Technical Data

	XXXX Bambino	XXS	Xs	S	M	L	XL
Seat width (cm)	8 - 16	12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)	11 x 6	14 x 6	18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)	220	220	220	220	280	280	280
Predetermined breaking point (at approx. kg)	40	40	60	60	80	80	80

Ratchet size XXXS - XXS	1,4 x 0,2 x 14,0 cm	<p>Product is latex-free</p> 	
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS	B-14-STAR-BL		
Buckle size XXS	B-14-GO-BL		
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			

PELVI.LOC® PL-3DSR Reverso

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis. On the moulded seat remain attached the ratchets. Remove the pads to implement.

Scope of application

For custom made seats, moulded seat.

Special features

The timing belts remain attached to the side of the seat pan. To transfer only remove the pads.

Fine adjustment of the pads to each other can be adjusted via the belt.

Article numbers

PL-3DSR-XXS for SW 12-20 cm

PL-3DSR-XS for SW 16-24 cm

PL-3DSR-S for SW 20-32 cm



PL-3DSR-M for SW 30-40 cm

PL-3DSR-L for SW 36-46 cm

PL-3DSR-XL for SW 44-52 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)		12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)		15 x 7	18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)		220	220	220	280	280	280
Predetermined breaking point (at approx. kg)		60	60	60	80	80	80

Ratchet size XXXS - XXS	1,9 x 0,2 x 40,0 cm	<p>Product is latex-free</p> 	
Ratchet size XS - S	1,9 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS	B-19-ERGO-BL		
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap XXS - XS - S	V-LR-19-BL-110 and V-LR-19-BL-140		
Hole-strap M - L - XL	V-LR-24-BL-110 and V-LR-24-BL-140		

PELVI.LOC® PL-3DSRR Reverso with ratchets

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Pelvis-retraction-belt with reversed ratchet

Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis.

The 3D formed shape is bolted firmly to the timing belt and strap hole.

Buckle belt **R-19/24-BL-400** and pressure-buckle **B-19/24-Star** are included.

Amendment

Hole-strap - set, available in a set of two

- Width: 19/24 mm
- Length: 115 mm

Hole-strap with buckle and clamping clip for mounting onto the frame. Also available in different sizes and diameters.

Example:

- Hole – strap 24 mm wide x 115 mm long and B-24-ERGO and clamping clip RS-25

Article numbers

LR-R-19-BL/B-19-ERGO-TR

LR-R-24-BL/B-24-ERGO-TR



Scope of application

For custom made seats and moulded seat.

Special features

The buckles are mounted on the side of the moulded seat or buggy. The pelvis-retraction-belt with the ratchet slides sideways into the buckles.

To improve positioning of the belt, adjustments can be made via the hole-strap.

Article numbers

PL-3DSRR-XS for width 16-24 cm



PL-3DSRR-S for width 20-32 cm

PL-3DSRR-M for width 30-40 cm

PL-3DSRR-L for width 36-46 cm

Technical Data

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)			16 - 24	20 - 32	30 - 40	36 - 46	
Padding size (cm)			18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	
Weight (g)			220	220	280	280	
Predetermined breaking point (at approx. kg)			60	60	80	80	

Ratchet size XXXS - XXS		<p>Product is latex-free</p> 	
Ratchet size XS - S	1,9 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-STAR-BL		
Buckle size M - L - XL	B-24-STAR-BL		
Hole-strap			

PELVI.LOC® PL-3DWP Waterproof

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Closed-pore special material, impermeable and water repellent! The pelvic pads are 3-dimensionally shaped for optimal fixation of the pelvis without pressure on the anterior superior spine. With easy-to-use stainless steel ratchet buckles, the fit of the strap is adjusted. The perforated belts on the frame can be adjusted up to 8cm in length and thus provide variable space for different situations. To secure the wheelchair user, the pelvic pads are connected in the middle via a toothed belt.



PL-3DWP-WH-L: 36-46 cm, weiß/white



PL-3DWP-BL-L: 36-46 cm, schwarz/black

Amendment

Pipe clamps for attaching the belt to the wheelchair frame available in the following diameters

Article numbers

RS-20 Tube diameter 20 mm
RS-22 Tube diameter 22 mm
RS-25 Tube diameter 25 mm
RS-30 Tube diameter 30 mm
RS-32 Tube diameter 32 mm



Scope of application

For pelvic positioning in the shower wheelchair, bathing chair, active wheelchair, rehab buggy, seat shells and much more.

Article numbers, black

PL-3DWP-BL-XXS for SW 12-20 cm
PL-3DWP-BL-XS for SW 16-24 cm
PL-3DWP-BL-S for SW 20-32 cm
PL-3DWP-BL-M for SW 30-40 cm
PL-3DWP-BL-L for SW 36-46 cm

Article numbers, white

PL-3DWP-WH-XXS for SW 12-20 cm
PL-3DWP-WH-XS for SW 16-24 cm
PL-3DWP-WH-S for SW 20-32 cm
PL-3DWP-WH-M for SW 30-40 cm
PL-3DWP-WH-L for SW 36-46 cm

Special feature

Facet-like thinning of the padding material in the transition from the pad to the fastening strap. Prevents lateral pressure points on the pelvis.

Ratchets, black



V-R-19-BL-ERGO-170
V-R-24-BL-ERGO-170

Ratchets, white

V-R-19-WH-ERGO-170
V-R-24-WH-ERGO-170

Technical Data

	XXXX Bambino	XXS	XS	S	M	L	XL
Seat width (cm)		12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)		14 x 6	18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)		220	220	220	280	280	280
Predetermined breaking point (at approx. kg)		40	60	60	80	80	80

Ratchet size XXXS - XS	1,9 x 0,2 x 17,0 cm	 <p>Product is latex-free</p>	 <p>LIZENZIERTER PARTNER</p> <p>REG. NR. DE3914274696739-V</p>
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			

PELVI.LOC® PL-ASY Asymmetrical Strap

PELVI.LOC® Seating systems (seat trousers)

HMV.Nr.: 18.99.99.0905



Version

Separately adjustable 3D-Contour pads take the Asymmetric sitting position on.

The resulting thrust - pressure of the pelvis forward is of the 3D - Contour taken the pads.

Special features

Prevents forward sliding of the basin, especially in reclining backrests and tilt.

Significant reduction in the pressure step by the 3D-Contour of the supports.

Easy handling due to pressure buckle **B-24-Star** star and Timing **R-24-BL-ERGO-400**.

Ratchet and buckle covers, prevents injuries.

Article numbers

PL-ASY-M Gr.M for width 36cm to 42cm

PL-ASY-S for SW 26-35 cm

PL-ASY-M for SW 36-39 cm


PL-ASY-L for SW 40-46 cm

Scope of application

For Asymmetric pelvic positioning in the seat, wheelchair or multifunction wheelchair.

Technical Data

	XXXS	XXS	XS	S	M	L	XL
Seat width (cm)				26 - 35	36 - 39	40 - 46	
Padding size (cm)				19 x 7,5	29 x 11	36 x 13,5	
Weight (g)				220	300	400	
Predetermined breaking point (at approx. kg)				60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	
Ratchet size XS - S	2,4 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XS - S	B-24-STAR-BL		
Buckle size M - L - XL	B-24-STAR-BL		

PELVI.LOC® PL-CSS Pelvi Contour Seating System

PELVI.LOC® Seating systems (seat trousers)

HMV.Nr.: 26.99.02.2000



Version

Pelvic retraction seating system.

Colour: black

The resulting thrust - pressure of the pelvis forward is of the 3D-Contour taken the pads. The padded seat prevents pressure.

Amendment

Easy and fast assembling the Zahriemens and webbing on the frame with a cable tie, belt loop or pipe clamp.



Scope of application

For pelvic positioning in a wheelchair or multifunction - Wheelchair.

Special features

Prevents forward sliding of the basin, especially in reclining backrests and tilt.

Significant reduction in the pressure step by the 3D-Contour of the supports.

Easy handling due to pressure buckle **B-24-Star** star and Timing **R-24-BL-ERGO-400**.

Ratchet and buckle covers, prevents injuries.

Article numbers



PL-CSS-S for SW 26-35 cm

PL-CSS-M for SW 36-39 cm

PL-CSS-L for SW 40-46 cm

Technical Data

	XXXS	XXS	XS	S	M	L	XL
Seat width (cm)				26 - 35	36 - 39	40 - 46	
Padding size (cm)				19 x 7,5	29 x 11	36 x 13,5	
Weight (g)				260	300	400	
Predetermined breaking point (at approx. kg)				60	80	80	

Ratchet size XXXS - XXS		<p>Product is latex-free</p> 	
Ratchet size XS - S	2,4 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS			
Buckle size XS - S	B-24-STAR-BL		
Buckle size M - L - XL	B-24-STAR-BL		

PELVI.LOC® PL-GF Genuframe

PELVI.LOC® Lower Leg Belts

HMV.Nr.: 18.99.99.0905



Version

The two-dimensionally shaped pads are set in front via the perforated belt. Timing belts **R24-BL-RNI-305** and buckle **B24-AE-BL** provide high load capacity, also available with other buckles and toothed belts.

The lower leg strap is individually adjusted to the sitting position by means of the ratchet buckle and can be adjusted to the nearest millimeter. Likewise, the strap can also be opened quickly via the ratchet buckle.

Amendment

With other buckles B24-star, B19-GO available.

Clamping clips for wheelchair frame.

Article numbers

RS-20 Tube diameter 20 mm

RS-22 Tube diameter 22 mm

RS-25 Tube diameter 25 mm

RS-30 Tube diameter 30 mm

RS-32 Tube diameter 32 mm



Scope of application

Pelvi.Loc - lower leg - belt prevents forward sliding of the pelvis in a wheelchair.

Special Features

Also with other buckles, e.g. (B24 Star, B19-GO) available.

Article numbers

PL-UG-XXS for SW 12-20 cm

PL-UG-XS for SW 16-24 cm

PL-UG-S for SW 20-32 cm



PL-UG-M for SW 30-40 cm

PL-UG-L for SW 36-46 cm

PL-UG-XL for SW 44-52 cm

Technical Data

	XXXX Bambino	XXS	XS	S	M	L	XL
Seat width (cm)		12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)		12 x 7	18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)		200	220	220	280	280	300
Predetermined breaking point (at approx. kg)		60	60	60	80	80	80

Ratchet size XXXS - XXS	1,9 x 0,2 x 40,0 cm	<p>Product is latex-free</p> 	
Ratchet size XS - S	1,9 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS	B-19-ERGO-BL		
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			

PELVI.LOC® PL-UG Lower Leg Belt

PELVI.LOC® Lower Leg Belts

HMV.Nr.: 18.99.99.0905



Version

The two-dimensionally shaped pads are set in front via the perforated belt. Timing belts **R24-BL-RNI-305** and buckle **B24-AE-BL** provide high load capacity, also available with other buckles and toothed belts.

The lower leg strap is individually adjusted to the sitting position by means of the ratchet buckle and can be adjusted to the nearest millimeter. Likewise, the strap can also be opened quickly via the ratchet buckle.

Amendment

With other buckles B24-star, B19-GO available.

Clamping clips for wheelchair frame.

Article numbers


RS-20 Tube diameter 20 mm

RS-22 Tube diameter 22 mm

RS-25 Tube diameter 25 mm

RS-30 Tube diameter 30 mm

RS-32 Tube diameter 32 mm



Scope of application

Pelvi.Loc - lower leg - belt prevents forward sliding of the pelvis in a wheelchair.

Special Features

Also with other buckles, e.g. (B24 Star, B19-GO) available.

Article numbers

PL-UG-XXXX for SW 8-16 cm

PL-UG-XXS for SW 12-20 cm

PL-UG-XS for SW 16-24 cm

PL-UG-S for SW 20-32 cm



PL-UG-M for SW 30-40 cm

PL-UG-L for SW 36-46 cm

PL-UG-XL for SW 44-52 cm

Technical Data

	XXXX Bambino	XXS	XS	S	M	L	XL
Seat width (cm)	8 - 16	12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)	11 x 6	14 x 6	18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)	220	220	220	220	280	280	280
Predetermined breaking point (at approx. kg)	40	40	60	60	80	80	80

Ratchet size XXXS - XXS	1,4 x 0,2 x 14,0 cm	<p>Product is latex-free</p> 	
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS	B-14-STAR-BL		
Buckle size XXS	B-14-GO-BL		
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		
Hole-strap			

PELVI.LOC® PL-UGR Lower Leg Belts reversed

PELVI.LOC® Lower Leg Belts

HMV.Nr.: 18.99.99.0905



Version

The two-dimensionally shaped pads are set in front via the perforated belt. Timing belts **R24-BL-RNI-305** and buckle **B24-AE-BL** provide high load capacity, also available with other buckles and toothed belts.

The lower leg strap is individually adjusted to the sitting position by means of the ratchet buckle and can be adjusted to the nearest millimeter. Likewise, the strap can also be opened quickly via the ratchet buckle.

Amendment

With other buckles B24-star, B19-GO available.

Clamping clips for wheelchair frame.

Article numbers


RS-20 Tube diameter 20 mm

RS-22 Tube diameter 22 mm

RS-25 Tube diameter 25 mm

RS-30 Tube diameter 30 mm

RS-32 Tube diameter 32 mm



Scope of application

Pelvi.Loc - lower leg - belt prevents forward sliding of the pelvis in a wheelchair.

Special Features

Also with other buckles, e.g. (B24 Star, B19-GO) available.

Article numbers

PL-UGR-XXS for SW 12-20 cm

PL-UGR-XS for SW 16-24 cm

PL-UGR-S for SW 20-32 cm

PL-UGR-M for SW 30-40 cm

PL-UGR-L for SW 36-46 cm

PL-UGR-XL for SW 44-52 cm

Technical Data

	XXXX Bambino	XXS	XS	S	M	L	XL
Seat width (cm)		12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)		14 x 6	18 x 7	21 x 8	25,5 x 8	26,5 x 8,5	28,5 x 9
Weight (g)		220	220	220	280	280	280
Predetermined breaking point (at approx. kg)		40	60	60	80	80	80

Ratchet size XXXS - XXS	1,9 x 0,2 x 40,0 cm	<p>Product is latex-free</p> 	
Ratchet size XS - S	1,9 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XXS	B-19-STAR-BL		
Buckle size XS - S	B-19-STAR-BL		
Buckle size M - L - XL	B-24-STAR-BL		
Hole-strap			

PELVI.LOC® PL-KG-GO Ankle Strap

PELVI.LOC® Foot Positioning



HMV.Nr.: 18.99.99.0305

Version

Serrated belt


Serrated belt **R-19-BL-VI-210 / 275 / 300** and Ratchet **B-19-GO-BL** incl. Accessory Kit WS-SET.

With timing belt cover after the buckle.

Included

WS-Set:

- Foot plate accessory kit M5 x 20 mm bolts with lock nuts, brackets and buckles.



Scope of application

Enables fast and easy location of foot /feet to chair. Specially designed to prevent injury to feet when not wearing shoes.

Article numbers

PL-KG-GO-S Circumference 17-23 cm



PL-KG-GO-M Circumference 20-26 cm

PL-KG-GO-L Circumference 24-30 cm

PL-KG-GO-XL Circumference 26-34 cm



Technical data

	S	M	L	XL	Product is latex free 	 LIZENZIERTER PARTNER <small>REG. NR. DE5914274696739-V</small>
Circumference (cm)	17 - 23	20 - 26	24 - 30	26-34		
Toothed strap	R-19-VI-BL-210	R-19-VI-BL-275	R-19-VI-BL-300	R-19-VI-BL-300		

PELVI.LOC® PL-RR Instep Strap

PELVI.LOC® Foot Positioning

HMV.Nr.: 18.99.99.0305

Foot straps

Padded instep strap for securing feet to wheelchair foot plate with 24 mm ratchet and serrated band.



Special features

Delivery always in pairs!
In many buckles and toothed belt variants possible!

Article numbers


- **PL-RR-XS**: Length min. 19 cm **B14-GO-BL** and **R14-GO-110**, Length max. 27 cm
- **PL-RR-S**: Length min. 20 cm **B24-LU-JU** and **R24-BL-170**, Length max. 27 cm
- **PL-RR-M**: Length min. 24 cm **B24-LU-JU** and **R24-BL-195**, Length max. 31 cm
- **PL-RR-L**: Length min. 24 cm **B24-LU-JU** and **R24-BL-225**, Length max. 36 cm
- **PL-RR-XL**: Length min. 29 cm **B24-LU-JU** and **R24-BL-225**, Length max. 51 cm
- **PL-RR-GO-S**: Length min. 20 cm **B19-GO-BL** and **R19-GO-BL-350**, Length max. 27 cm
- **PL-RR-GO-M**: Length min. 24 cm **B19-GO-BL** and **R19-GO-BL-350**, Length max. 31 cm

Included

BS-BL-25-SET

Mounting Bracket Set with screws and nuts.

Always in pairs!



Scope of application

Padded instep strap for securing feet to wheelchair foot plate with 24 mm ratchet and serrated band.
Allows for quick release and securing with one hand.



Version PL-RR

Upholstery size above everything:

- **Size: XS-S** 17,5 cm x 5,0 cm, overall length 21,0 cm
- **Size: M** 17,5 cm x 5 cm, overall length 23,5 cm
- **Size: L** 17,5 cm x 5 cm, overall length 23,5 cm, same cushion size wie Size M



Version PL-RR

Upholstery size above everything:

- **Size: XL** 24,0 cm x 7,0 cm, overall length 35,0 cm



Version PL-RR-GO

Upholstery size above everything:

- **Size: S** 20,0 cm, overall length 27,0 cm
- **Size: M** 24,0 cm, overall length 31,0 cm

Technical data

PL-RR	XS	S	M	L	XL
Padding size (cm)	17,5 x 5	17,5 x 5	25,5 x 8	25,5 x 8	24 x 7
Total length (cm) minimal	19	20	24	24	29
Total length (cm) maximal	27	27	31	36	51

Product is latex free





Technik die verbindet
Technology that connects

Toothed strap	R-14-GO-110	R-24-BL-170	R-24-BL-195	R-24-BL-225	R-24-BL-225
Buckle	B-14-GO-BL	B-24-LU-JU	B-24-LU-JU	B-24-LU-JU	B-24-LU-JU

PL-RR-GO	S	M		
Padding size (cm)	17,5 x 5	25,5 x 8		
Total length (cm) minimal	20	24		
Total length (cm) maximal	27	31		
Toothed strap	R-19-GO-BL-350	R-19-GO-BL-350		
Buckle	B-19-GO-BL	B-19-GO-BL		

Product is latex free

PELVI.LOC® / Torso.Loc Two in One TL-3DSR & TL-W

PELVI.LOC® / Torso.Loc Upper Body Positioning

HMV.Nr.: 18.99.99.0906



Version

The Torso.Loc offers an exclusive combination of upper body positioning and pelvic retraction harness. The neoprene vest has a high restoring force and is individually adjustable in the shoulder area and sternum for a secure positioning of the upper body. Vest and belt are connected via a soft and tension-resistant belt in the upper abdomen area. The 3D Pelvic Retraction Belt provides optimal pelvic fixation without pressure on the anterior superior spine. The Torso.Loc ensures optimum positioning of the upper body and pelvis in the wheelchair.



Special features

Two in one: The Torso.Loc system is characterized by simple and quick handling: the combination of the neoprene upper body vest and the TL-3DSR pelvic retraction belt reduces several work steps to one step for optimal positioning.

Article numbers: TL-3DSR-XS / TL-W-S

Scope of application

For upper body and pelvic positioning in the wheelchair, the seat shell or rehab buggy. When using the Torso.Loc system, with the PEG gavage lying, one must pay attention to a possible pressure load on the puncture site and the PEG tube.

Consisting of Pelvic retraction belt




Article numbers

TL-3DSR-XS for SB 16-24 cm
TL-3DSR-S for SB 20-32 cm
TL-3DSR-M for SB 30-40 cm
TL-3DSR-L for SB 36-46 cm

Consisting of neoprene vest





Pictogram neoprene vest



Article numbers	A	B	C	D	E	F
TL-W-S	31	22	12	16	4	12
TL-W-M	38	29	16	19	5	13
TL-W-L	47	37	19	22	6	14
TL-W-XL	55	45	24	26	6,5	17

Technische Daten TL-3DSR

	XS	S	M	L	Product is latex-free   LIZENZIERTER PARTNER REG. NR. DE5914274696739-V
Seat width (cm)	16 - 24	20 - 32	30 - 40	36 - 46	
Padding size (cm)	18 x 7	20 x 11	23 x 13	25,5 x 14	
Weight (g)	220	220	280	280	
Predetermined breaking point (at approx. kg)	60	60	80	80	
Ratchet	1,9 x 0,2 x 40,0 cm	1,9 x 0,2 x 40,0 cm	2,4 x 0,2 x 40,0 cm	2,4 x 0,2 x 40,0 cm	

PELVI.LOC® PL-OPS Upper Body Positioning made to measure

PELVI.LOC® / Torso.Loc Upper Body Positioning

HMV.Nr.: 18.99.99.0906



Version

Our upper body positioning system is individually adaptable. According to your measurements, we manufacture the system for safe sitting with upper body stabilization. For the car transport of children, adolescents and adults, the system offers simple and quickly applicable handling, and the upper body positioning can be installed in just a few simple steps.

OPS made to measure

assembled in a few easy steps!

Amendment

Safety buckle protects against unintentional opening of the body belt.



Article number: SF-2

Segufix magnetic closure with 2 magnetic closures incl. 2 security keys to open the system. Only on request and agreement possible!

Article number:

KSS-S-50 3 finger buckle

KSS-50 2 finger buckle

SKS-50-SF Lock with safety function (pin opening)

Scope of application

For the upper body - positioning in a wheelchair, powered wheelchair, car Seat.

Special features

Upper body - positioning - system to your specifications

- Body part with lateral stability
- Seat - belt individually adjustable in length by Velcro fastener
- Wide, padded shoulder straps
- Special – Hook/Loop fasteners with very high tensile force and Life

Article number: PL-OPS-Dimension



Article number: KSS-S-50

To open lock with three fingers - Nylon - Pressure



Article number: KSS-50

To open lock with two fingers - Nylon - Pressure



Article number: SKS-50-SF

Metal lock 50mm wide with safety function (Pin opening)

PELVI.LOC® PL-W

PELVI.LOC® / Torso.Loc Upper Body Positioning

HMV.Nr.: 18.99.99.0906



Version


Neoprene harness Classic without zip

Our neoprene harnesses are individually adjustable. For secure positioning of the upper body in the wheelchair.

Scope of application

For upper body positioning, stabilisation in wheelchairs, care wheelchairs.

Pictogram PL-W



Article numbers	A	B	C	D	E	F
PL-W-S	31	22	12	16	4	12
PL-W-M	38	29	16	19	5	13
PL-W-L	47	37	19	22	6	14
PL-W-XL	55	45	24	26	6,5	17

Special features

Material: Neoprene 7 mm laminated, ensures a high restoring force in the shoulder area.

PELVI.LOC® PL-WR

PELVI.LOC® / Torso.Loc Upper Body Positioning

HMV.Nr.: 18.99.99.0906



Version


Neoprene harness Classic with zip

Our neoprene waistcoats are individually adjustable. For secure positioning of the upper body in the wheelchair.

Scope of application

For upper body positioning, stabilisation in wheelchairs, care wheelchairs.

Pictogram PL-WR



Article numbers	A	B	C	D	E	F
PL-WR-S	31	22	12	16	4	12
PL-WR-M	38	29	16	19	5	13
PL-WR-L	47	37	19	22	6	14
PL-WR-XL	55	45	24	26	6,5	17

Special features

Material: Neoprene 7 mm laminated, ensures a high restoring force in the shoulder area.