# PELVI.LOC® PL-2DA Active

PELVI.LOC® PL-2D Thigh Belts

HMV.Nr.: 18.99.99.0905



### Version

2D-belt for enhanced sports activities with enlarged pads - surface for improved pressure distribution.

The padded covers give maximum protection of participants and prevent buckle from releasing in use. Hole strap on frame adjustable up to 6 cm.

## **Special features**

Monoski racer with Martin Braxenthaler:

• PL-2DA-M as a chest strap and thigh strap



## **Amendment**

Tube clamps for attaching the belt to the wheelchair frame available in the following diameters:

#### Article numbers

RS-20: Tube diameter 20 mm RS-22: Tube diameter 22 mm RS-25: Tube diameter 25 mm RS-30: Tube diameter 30 mm



## Scope of application

- Safe driving with 2D pelvic positioning in a wheelchair.
- Can be used for thorax positioning in a standing device.
- Pelvic positioning with pressure distribution and stability in sports wheelchairs.

## Special features

#### Article numbers

PL-2DA-S for SW 20-32 cm PL-2DA-M for SW 30-40 cm PL-2DA-L for SW 36-46 cm

## **Technical Data**

	XXXS Bambino	XXS	XS	S	M	L	XL
Seat width (cm)				20 - 32	30 - 40	36 - 46	
Padding size (cm)				19 x 11	23 x 12	28 x 14	
Weight (g)				220	280	280	
Predetermined breaking point (at approx. kg)				60	80	80	

Ratchet size XXXS - XXS	
Ratchet size XS - S	2,4 x 0,2 x 17,0 cm
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm
Buckle size XXXS	
Buckle size XXS	
Buckle size XS - S	B-24-FLAT-LO-BL
Buckle size M - L - XL	B-24-ERGO-BL
Hole-strap	LR-BL-320

Produkt ist latexfrei



