

PELVI.LOC® PL-2DB

PELVI.LOC® PL-2D Thigh Belts

HMV.Nr.: 18.99.99.0905



Version

2DB belt suitable for thigh. The buckles are easily accessible

The ratchet will run in addition to a strap to prevent injuries. Hole strap on frame adjustable up to 6 cm. For fixing on to the frame use compatible clamping clips.

Performance

Cross country skis with:

- **PL-3DSRR-M** 3D moulded seat Reverso with twisted toothed belt and belt set
- **PL-2DB-XS** Thigh Belt Front



Amendment

Tube clamps for attaching the belt to the wheelchair frame available in the following diameters:

Article numbers

RS-20: Tube diameter 20 mm
RS-22: Tube diameter 22 mm
RS-25: Tube diameter 25 mm
RS-30: Tube diameter 30 mm
RS-32: Tube diameter 32 mm



Scope of application

- For pelvis positioning in a wheelchair and sports wheelchair

Article numbers

PL-2DB-XXS for SW 12-20 cm
PL-2DB-XS for SW 16-24 cm
PL-2DB-S for SW 20-32 cm
PL-2DB-M for SW SB 30-40 cm
PL-2DB-L for SW SB 36-46 cm
PL-2DB-XL for SW 44-52 cm



Technical Data

	XXXS	XXS	XS	S	M	L	XL
Seat width (cm)		12 - 20	16 - 24	20 - 32	30 - 40	36 - 46	44 - 52
Padding size (cm)		12 x 7	17 x 7	20 x 7	22 x 8	25,5 x 8	28,5 x 9
Weight (g)		200	220	220	280	280	280
Predetermined breaking point (at approx. kg)		40	60	60	80	80	8

Ratchet size XXXS - XXS	1,4 x 0,2 x 14,0 cm	Product is latex-free 	 LIZENZIIERTER PARTNER REG. NR. DE 591427/4696739-V
Ratchet size XS - S	1,9 x 0,2 x 17,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 17,0 cm		
Buckle size XXXS	B-14-STAR-BL		
Buckle size XXS	B-14-GO-BL		
Buckle size XS - S	B-19-ERGO-BL		
Buckle size M - L - XL	B-24-ERGO-BL		