

PELVI.LOC® PL-3DASR Activ reversed for moulded Seat

PELVI.LOC® PL-3D Pelvis Belts

HMV.Nr.: 26.99.02.2000



Version

Reversed pelvis retraction belt for moulded seat

Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis. 3-Dimensional shape gives maximum support of the pelvis without pressure on the spine "anterior superior".

Formed shape with ultra thin end at the moulded seat fixing point to avoid lateral pressure marks at the pelvis. The covered buckles help protect from injury and accidental opening during use.

Amendment

Pipe clamps for attaching the belt to the wheelchair frame available in the following diameters.

Article numbers

RS-20 Tube diameter 20 mm
RS-22 Tube diameter 22 mm
RS-25 Tube diameter 25 mm
RS-30 Tube diameter 30 mm
RS-32 Tube diameter 32 mm



Scope of application

For custom-made seating systems, moulded seat.

Special features



On the moulded seat remain attached the ratchets. Remove to improve positioning of the pelvis - retraction of the seat, with high spasticity.

Article numbers

PL-3DASR-XS for width 16-24 cm
PL-3DASR-S for width 20-32 cm
PL-3DASR-M for width 30-40 cm
PL-3DASR-L for width 36-46 cm

Technical Data

| | XXXS Bambino | XXS | XS | S | M | L | XL |
|--|--------------|-----|----------|---------|---------|-----------|----|
| Seat width (cm) | | | 16 - 24 | 20 - 32 | 30 - 40 | 36 - 46 | |
| Padding size (cm) | | | 17,5 x 9 | 20 x 11 | 23 x 13 | 25,5 x 14 | |
| Weight (g) | | | 220 | 220 | 280 | 280 | |
| Predetermined breaking point (at approx. kg) | | | 60 | 60 | 80 | 80 | |

| | | | |
|-------------------------|---------------------|---|--|
| Ratchet size XXXS - XXS | | Product is latex-free  |  LICENZIIERTER PARTNER REG. NR. DE5914274656739-V |
| Ratchet size XS - S | 1,9 x 0,2 x 40,0 cm | | |
| Ratchet size M - L - XL | 2,4 x 0,2 x 40,0 cm | | |
| Buckle size XXXS | | | |
| Buckle size XXS | | | |
| Buckle size XS - S | B-19-ERGO-BL | | |
| Buckle size M - L - XL | B-24-ERGO-BL | | |
| Hole-strap | | | |