

PELVI.LOC® PL-ASY Asymmetrical Strap

PELVI.LOC® Seating systems (seat trousers)

HMV.Nr.: 18.99.99.0905



Version

Separately adjustable 3D-Contour pads take the Asymmetric sitting position on.

The resulting thrust - pressure of the pelvis forward is of the 3D - Contour taken the pads.

Special features

Prevents forward sliding of the basin, especially in reclining backrests and tilt.

Significant reduction in the pressure step by the 3D-Contour of the supports.

Easy handling due to pressure buckle **B-24-Star** star and Timing **R-24-BL-ERGO-400**.

Ratchet and buckle covers, prevents injuries.

Article numbers

PL-ASY-M Gr.M for width 36cm to 42cm

PL-ASY-S for SW 26-35 cm

PL-ASY-M for SW 36-39 cm

PL-ASY-L for SW 40-46 cm

Scope of application

For Asymmetric pelvic positioning in the seat, wheelchair or multifunction wheelchair.

Technical Data

	XXXS	XXS	XS	S	M	L	XL
Seat width (cm)				26 - 35	36 - 39	40 - 46	
Padding size (cm)				19 x 7,5	29 x 11	36 x 13,5	
Weight (g)				220	300	400	
Predetermined breaking point (at approx. kg)				60	80	80	

Ratchet size XXXS - XXS		Product is latex-free 	
Ratchet size XS - S	2,4 x 0,2 x 40,0 cm		
Ratchet size M - L - XL	2,4 x 0,2 x 40,0 cm		
Buckle size XXXS			
Buckle size XS - S	B-24-STAR-BL		
Buckle size M - L - XL	B-24-STAR-BL		