

PELVI.LOC® PL-CSS Pelvi Contour Seating System

PELVI.LOC® Seating systems (seat trousers)

HMV.Nr.: 26.99.02.2000



Version

Pelvic retraction seating system.

Colour: black

The resulting thrust - pressure of the pelvis forward is of the 3D-Contour taken the pads. The padded seat prevents pressure.

Amendment

Easy and fast assembling the Zahriemens and webbing on the frame with a cable tie, belt loop or pipe clamp.







Scope of application

For pelvic positioning in a wheelchair or multifunction - Wheelchair.

Special features

Prevents forward sliding of the basin, especially in reclining backrests and tilt.

Significant reduction in the pressure step by the 3D-Contour of the supports.

Easy handling due to pressure buckle **B-24-Star** star and Timing **R-24-BL-ERGO-400**.

Ratchet and buckle covers, prevents injuries.

Article numbers

PL-CSS-S for SW 26-35 cm **PL-CSS-M** for SW 36-39 cm **PL-CSS-L** for SW 40-46 cm

Technical Data

| | XXXS | XXS | XS | S | M | L | XL |
|--|------|-----|----|----------|---------|-----------|----|
| Seat width (cm) | | | | 26 - 35 | 36 - 39 | 40 - 46 | |
| Padding size (cm) | | | | 19 x 7,5 | 29 x 11 | 36 x 13,5 | |
| Weight (g) | | | | 260 | 300 | 400 | |
| Predetermined breaking point (at approx. kg) | | | | 60 | 80 | 80 | |

| Ratchet size XXXS - XXS | |
|-------------------------|---------------------|
| Ratchet size XS - S | 2,4 x 0,2 x 40,0 cm |
| Ratchet size M - L - XL | 2,4 x 0,2 x 40,0 cm |
| Buckle size XXXS | |
| Buckle size XXS | |
| Buckle size XS - S | B-24-STAR-BL |
| Buckle size M – L – XL | B-24-STAR-BL |

Product is latex-free



